



## BD Accuri C6 Template

Detector Array	Position of PMT	Lambda of Band-Pass Filter	Fluorochrome (Example)
488 Laser (20mW)	FL1	530/30	FITC / GFP/ Alexa Fluor 488 / YFP / CFSE
	FL2	585/40	PE / PI/PE-Texas RED
	FL3	LP 670	Pe Cy5 / 7AAD / Pe Cy5.5 / PerCP / PerCP Cy 5.5 / PE Cy7
640 Laser (12.5mW)	FL4	675/25	APC / Alexa Fluor 647