

In collaboration with



WEBINARDate 24/03 at 11.00 AMDEMO ON-SITEDate 4/04 - 5/04BOOK YOUR SESSION

W8 PHYSICAL CYTOMETER

The W8 Physical Cytometer: Mass-Density and Weight as Novel Biophysical Parameters to Analyse and Sort 3D Cell Culture

flow-based technology, The which is gathering more and more attention since the official launch in 2021. combines sedimentation theory with customized image-based tracking analysis and software algorithms, to determine samples' mass density, weight and size.



SPEAKERS





Domenico Andrea Cristaldi PhD R&D Specialist and Tech Marketing at CellDynamics

Spartaco Santi PhD Head of Digital Microscopy Center at IOR